

Product Features

- Uses multiple 20-slide baskets
- Intuitive software makes it possible to perform separate protocols simultaneously with multiple baskets in process
- Units available with or without five independently heated positions
- 6 water stations, 26 reagent stations, 2 loading and unloading ports
- IVDD compliant

Powerful Software

- Allocates reagents, optimises the reagent plan and calculates the most efficient route for each protocol
- Graphically displays the progress of each basket along with a text summary on the unit's VGA monitor
- Comes complete with suggested staining protocols, which may be easily modified to the laboratory's procedures
- Stores as many as 50 staining procedures in memory

Unique Design

- Two level design saves bench space
- Lower-level staining stations are on a turntable for easy access to all dishes from the front of the machine

Specifications

Model

Dimensions H x W x D mm

Weight kg (lbs)

Safety

- Continues even if power fails
- Handles fume control with easy-to-change carbon filters or through connection to an optional external vent
- Internal flow regulators and anti-siphon devices on the water supply lines protect the water supply
- Smart software keeps the unit's robotic arm from entering user access areas

Shandon Varistain Gemini

Stainer	800 x 720 x 710
Power Supply	160 x 350 x 468
Stainer	80 (176)
Power Supply	10 (22)

Ordering Information

Cat No.

A78000002

A78000012

Description

Shandon Varistain Gemini automatic slide stainer without heat, includes monitor
 Shandon Varistain Gemini automatic slide stainer with heat, includes monitor

Accessories

Cat No.

Included with system

A78010201

A78010466

A78010487

9990610

AP14185

AP14160

A78030174

Optional Accessories

A78010403

A78010404

Description

3-well water stations x 2
 Staining basket/carriers 2 x 5 pk
 Spare staining troughs 3 pk
 Activated charcoal filter
 Monitor cable
 Water inlet hose
 Screened data cable

Consul transfer clip tool
 Sakura basket carrier

