

The Ambassador 100 / 125 / 158



Features

- * Available in 100, 125 or 158 litre capacity formats.
- * MP 24 controller (See page 34)
- * Robust two bolt door closure system.
- * Proven reliable design.
- * Total service accessibility.
- * Integral drip tray
- * Total after sales service.
- * CE marked.
- * Ambassador 100 is single phase electrically heated.
- * Ambassador 125/158 available in electrically heated versions or from a piped steam source.
- * Electrically heated versions generate steam from a water well below the shelf.

Optional Accessories

- * Loading baskets.
- * Waste discard system.
- * High-speed fan assisted cooling.
- * Drain line condenser for connection to plastic drains.
- * Full range of printers and cycle recorders.
- * Condense collector - eliminates the need for a drain.
- * Vacuum air removal system.
- * Automatic water fill.
- * Separate steam generator mounted below the chamber.
- * Pre-heated steam generator.
- * Loading tray.

Specifications

Part No.	C/A/010100Amb	C/A/010125Amb	C/A/010158Amb
Chamber Capacity litres	100	125	158
Vessel Dims (mm)	460Ø x 610dp	460Ø x 760dp	460Ø x 960dp
Overall Dims (mm) wxdxh	650x1160x1240	650x1160x1240	650x1310x1240
Power Requirements (Elec)	6kW, sing phase	9kW, 3 phase	9kW, 3 phase
Power Requirements (Stm)	N/A	Single phase	Single phase
Cooling	-----Fan Assisted Cooling Optional-----		
Drain Line	-----Preferred > 22mm-----		
Water Supply	-----Optional 15mm-----		