



Digital controlled range of water baths.
 Double walled metallic construction, stainless steel (18/10) interior, exterior wall painted in epoxy and treated with antioxidant.
 Heating by means of element shielded in stainless steel AISI 316.
 Temperature regulation through digital thermostat.
 Equipped with safety thermostat with manual re-set according to DIN 12877 class 2.
 Control Panel with:
 -Illuminated start-up switch.
 -Pilot light indicating electric resistance operation.
 -Pilot light indicating safety thermostat operation.
 -Temperature adjustment.
 Manufactured according to EC directives.

BD Series

Maximum temperature - 100°C

Part No.	Capacity (litres)	Measurements W x H X D (cm)	Ext. Measurements W x H X D (cm)	Power (w)	Probe type	Probe Precision
C/B/010034	5	29 x 15 x 11	37 x 20 x 28	500	PTC	±2%
C/B/010035	12	32 x 15 x 29	39 x 20 x 11	1200	PTC	±2%
C/B/010036	20	49 x 15 x 29	56 x 20 x 11	1700	PTC	±2%

Maximum temperature - 200°C

Part No.	Capacity (litres)	Measurements W x H X D (cm)	Ext. Measurements W x H X D (cm)	Power (w)	Probe type	Probe Precision
C/B/010037	5	29 x 15 x 14	37 x 20 x 28	700	pt100	±0.5%
C/B/010038	12	32 x 15 x 29	39 x 20 x 44	1500	pt100	±0.5%
C/B/010039	20	49 x 15 x 29	56 x 20 x 44	2000	pt100	±0.5%

Accessories

Drainage tap; **Part No. C/B/010268**

Racks for 98 test tubes of 10-12 mm diameter; **Part No. C/B/010272**

Racks for 60 test tubes of 14-16 mm diameter; **Part No. C/B/010273**

Racks for 50 test tubes of 18-20 mm diameter; **Part No. C/B/010274**

Level liquid indicator; **Part No. C/B/010271**

For Baths with Part No.	Lid with holes and reduction disks 40, 60, 80 (mm)	Top Cover	Blank Cover	Element Tray cover
C/B/010034, C/B/010037	C/B/010332 two places	C/B/010275	C/B/010278	C/B/010281
C/B/010035, C/B/010038	C/B/010333 four places	C/B/010276	C/B/010279	C/B/010282
C/B/010036, C/B/010039	C/B/010334 six places	C/B/010277	C/B/010280	C/B/010283