

Steddy Range with Different Eyepieces and Stands

C/M/010043	Optical head stereoscopic Bino STEDDY-B, with eyepieces 10x
C/M/010044	Optical head stereoscopic Trino STEDDY-T, with eyepieces 10x
C/M/010045	Focusing device to adapt bino-/trinocular head to universal stand incl. vertical column 25mm
C/M/010046	Universal swivel stand + table clamp
C/M/010047	Universal swivel stand + rectangular base
C/M/010048	Universal stand + round base

Accessories

C/M/010148	Cable with US plug and Euro connector, black 2m
C/M/010149	Cable with UK plug and Euro connector, black 2m
C/M/010150	Cleaning set: lens cleaning tissues (20 pcs) + solution (75ml) + brush
C/M/010180	Eyepiece micrometer single 10x/18mm with cross line scale 0-14mm : 0,1mm
C/M/010181	Eyepiece pair 5x/22mm wide-field
C/M/010182	Eyepiece pair 15x/13mm wide-field
C/M/010183	Eyepiece pair 20x/10mm wide-field
C/M/010184	Additional lens 0,5x negative lens: increases working distance and increases field of view
C/M/010185	Additional lens 0,75x negative lens: increases working distance and increases field of view
C/M/010186	Additional lens 1,5x positive lens: decreases working distance and decreases field of view
C/M/010187	Condenser darkfield
C/M/010188	Jewelry clamp
C/M/010189	Polarising set
C/M/010190	Microscope adapter for single lens reflex camera including projection lens 2x
C/M/010191	Microscope adapter for 1/3 and 1/2 inch CCD-camera , including C-mount
C/M/010430	Microscope adapter ring illuminator for Steddy

Spare Parts

C/M/010220	Fuse 5 x 20/500 MA (10 pcs)
C/M/010278	Eyepiece pair 10x/20mm wide-field
C/M/010279	Filter blue Ø 35mm
C/M/010282	Eyecup pair, rubber
C/M/010283	Stage plate, frosted ø 80mm
C/M/010284	Stage plate, black/white ø 80mm
C/M/010285	Stage clamp, pair
C/M/010286	Projection lens 2x for photomicrography
C/M/010367	Halogen reflector bulb 12V/10W (10 pcs) for episcopic light
C/M/010368	Bulb halogen 12V/10W/G4 (10 pcs) for diasopic light

