

Twister Range with Different Eyepieces and Stands

C/M/010051	Microscope Stereo TWISTER Zoom with Flat base no illumination
C/M/010052	Optical head stereoscopic Bino TWISTER Zoom, with eyepieces 10x
C/M/010053	Optical head stereoscopic Bino TWISTER Step, with eyepieces 10x
C/M/010054	Focusing device Twister to adapt binocular head to universal stand incl. vertical column 30mm
C/M/010055	Universal swivel stand + table clamp
C/M/010056	Universal swivel stand + rectangular base
C/M/010057	Universal stand + round base

Accessories

C/M/010148	Cable with US plug and Euro connector, black 2m
C/M/010149	Cable with UK plug and Euro connector, black 2m
C/M/010150	Cleaning set: lens cleaning tissues (20 pcs) + solution (75ml) + brush
C/M/010431	Eyepiece micrometer single 10x/18mm with cross line scale 0-14mm : 0,1mm
C/M/010432	Eyepiece pair 5x/22mm wide-field
C/M/010434	Eyepiece pair 15x/16mm wide-field
C/M/010435	Eyepiece pair 20x/13mm wide-field
C/M/010436	Eyepiece pair 25x/12mm wide-field
C/M/010437	Additional lens 0,5x negative lens: increases working distance and increases field of view
C/M/010438	Additional lens 0,75x negative lens: increases working distance and increases field of view
C/M/010439	Additional lens 1,5x positive lens: decreases working distance and decreases field of view
C/M/010444	Additional lens 2,0x positive lens: decreases working distance and decreases field of view
C/M/010440	Microscope adapter for SLR camera including projection lens
C/M/010441	Microscope adapter for 1/3 inch CCD-camera with 0,3x lens, including C-mount
C/M/010442	Microscope adapter for 1/2 inch CCD-camera with 0,5x lens, including C-mount
C/M/010443	Beam splitter, Trinocular tube

Spare Parts

C/M/010220	Fuse 5 x 20/500 MA (10 pcs)
C/M/010369	Filter blue Ø 40mm
C/M/010370	Bulb halogen reflector 6v/12w for episcopic light
C/M/010371	Bulb fluorescent tubular 6v/15w diascopic light
C/M/010372	Eyecup pair, rubber
C/M/010373	Stage plate, frosted Ø 95mm
C/M/010374	Stage plate, black/white Ø 95mm
C/M/010375	Stage clamp, pair
C/M/010376	Eyepiece pair 10x/23mm super wide-field for Twister S
C/M/010433	Eyepiece pair 10x/22mm super wide-field for Twister Z