

For Quality Results

- Temperature range up to 220°C
- Analog controller and dial temperature gauge for ease of use
- Robust construction of epoxy coated high quality steel and aluminum
- Triple insulated body ensures excellent chamber homogeneity and minimises power consumption
- Wide range of options available including stainless steel interior, digital display, Class 3:1 safety thermostat and timer
- Two large capacity models available
- Reliable performance with electronic regulation
- Ideal for glassware sterilisation of large loads
- Delivered with two gridded shelves



Specifications

Model

	EU500	EU720
Capacity ltrs (cu. ft)	493 (17.3)	721 (25.3)
No. Shelves Std/Max	2/19	2/28
Temperature Range	60°C to 220°C	60°C to 220°C
Exterior Dimensions (H x W x D) mm (in)	1385 x 1260 x 630 (54.6 x 49.7 x 24.8)	1800 x 1260 x 630 (71 x 49.7 x 24.8)
Interior Dimensions (H x W x D) mm (in)	905 x 1100 x 495 (54 x 49.5 x 24.7)	1325 x 1100 x 495 (52 x 43 x 19.4)
Net Weight kg (lbs)	195 (432)	215 (477)

Ordering Information

Cat No.

46151002

46151004

46161002

46161004

Description

EU500 Large Capacity Oven (air extraction)
 EU500 Large Capacity Oven (forced ventilation)
 EU500 Large Capacity Oven (air extraction)
 EU500 Large Capacity Oven (forced ventilation)

Power Requirements

230V / 50Hz
 230V / 50Hz
 230V / 50Hz
 230V / 50Hz

Jouan EU Series Options and Accessories

Options and Accessories

Cat No.

Options must be ordered with oven

00002000

00050001

00060001

00000200

00000400

00000800

00010010

Description

Digital display - for EU500-720
 Stainless steel interior - for EB/EU500
 Stainless steel interior - for EB/EU720
 12 hour timer*
 Delayed start timer*- for EB500-950/EU500-720
 Dual temperature timer* - for EB500-950/EU500-720
 Cable access port - for EB500-950/EU500-720

Accessories

41390004

41390007

41390005

Stainless steel gridded shelf - for EU500-720
 Wide stainless steel shelf - for EU500-720
 Humidity tray - for EB500-950/EU500-720

* Only one timer at a time