

Designed for a long trouble-free life, the Jouan EB Series large capacity incubators are supremely reliable. Their conception with powerful heaters and electronic regulation makes them reliable with temperature overshoot and fluctuations reduced to a minimum.

Product Features

- Easy to use with temperature range up to 70°C, simple analog controller and dial temperature gauge
- Robust construction and triple insulated body that ensure good chamber homogeneity and minimises power consumption
- Versatility is ensured by multiple options to adapt to specific applications
- Selection of three capacities to fit all large loads and space constraints



Specifications

Cat No.	41353000	41363000	41373000
Capacity cu ft (ltrs)	17.3 (493)	25.3 (721)	33.2 (947)
Temperature Regulation Range	Ambient +5°C to 70°C	Ambient +5°C to 70°C	Ambient +5°C to 70°C
Exterior Dimension mm	1385 x 1260 x 630	1800 x 1260 x 630	2220 x 1260 x 630
Interior Dimension mm	905 x 1100 x 495	1325 x 1100 x 495	1740 x 1100 x 495
Net Weight kg (lbs)	195 (432)	215 (477)	235 (521)

Ordering Information

Cat No.	Description	Power Requirements
41353000	EB500 Large Capacity Bacteriological Incubator	230 V / 50 Hz
41363000	EB720 Large Capacity Bacteriological Incubator	230 V / 50 Hz
41373000	EB950 Large Capacity Bacteriological Incubator	230 V / 50 Hz

Options and Accessories

Options must be ordered with the Incubator

00100000	Digital display for EB500-950
00100004	Forced ventilation for EB500-950
00050001	Stainless steel interior - for EB/EU500
00060001	Stainless steel interior - for EB/EU720
00070001	Stainless steel interior for EB950
00000200	12 hour timer*
00000400	Delayed start timer*- for EB500-950/EU500-720
00000800	Dual temperature timer* - for EB500-950/EU500-720
00600020	Quadruple internal sockets for EB500-950
00010010	Cable access port - for EB500-950/EU500-720

Accessories

41390001	Perforated galvanised steel shelf + epoxy coated for EB500-950
41390002	Perforated galvanised steel shelf for EB500-950
41390003	Lowenstein shelf (63 tubes) for EB500-950
41390010	Lowenstein shelf with 3 drawers (3 x 21 tubes) for EB500-950
41390009	Lowenstein shelf for 170 mm tubes (2 x 21 tubes) for EB500-950
41390005	Humidity tray - for EB500-950/EU500-720
41390024	7 day battery operated temperature chart recorder 0-75°C (delivered with 100 recorder discs)
41390025	100 recorder discs

Note: * Only one timer at a time