

The Jouan Artic cooled incubators are designed to meet the needs of incubation below and around ambient temperature. These incubators are very easy to use with all controls and displays on the front panel. They can maintain very accurate temperatures with good temperature homogeneity guaranteed by the internal circulation fan and well insulated chamber.

### Product Features

- Ideal for all applications below and around ambient temperatures as well as standard incubation temperatures up to 40°C
- Versatility of use in precise environment: they are designed to accurately hold temperatures at and below ambient. The critical temperatures of 4°C and 20°C are precisely maintained. In addition, a temperature range of up to 40°C allows use as a standard incubator.
- Two capacities: 130L (Mini Artic) and 260L (Maxi Artic)
- Accurate temperature control over +2°C to + 40°C with Pt probe
- Large storage surface with a maximum of 15 or 35 shelves
- Optimised use of space; reversible door
- Excellent homogeneity with forced ventilation.



### Applications

During summer, maintaining a room temperature can be a problem. The

Jouan Artic Cooled incubators are ideal for holding 20°C all year round. The Jouan Artic cooled incubators meet the

needs of BOD, tissue culture, serum storage, water analysis, cooled blooded animal tests, fungal analysis.

### Specifications

#### Model

Capacity ltrs (cu ft)

Temperature Regulation Range

Exterior Dimensions mm

Interior Dimensions mm

Net Weight kg (lbs)

#### Mini Artic

130 (4.5)

2°C to 40°C

950 x 602 x 600

702 x 513 x 433

45 (99.7)

#### Maxi Artic

260 (9.1)

2°C to 40°C

1673 x 602 x 600

1418 x 513 x 433

64 (141.8)

### Ordering Information

#### Cat No.

11326301

11326351

#### Description

Mini Artic Cooled Incubator

Maxi Artic Cooled Incubator

#### Power Requirements

230 V / 50 Hz

230 V / 50 Hz

### Options and Accessories

Options must be ordered with the Incubator

0000011

0000010

Internal Power Socket

Access Port

### Accessories

86000872

Extra Shelf

Note:- The Jouan Cooled Incubators are delivered with the following standard accessories : 3 x shelves (Mini), 5 x shelves (Maxi)

