



Buffer Solutions

Orion present buffers in different formats according to customer choice. Most popular are the twin necked bottles which allow calibration to be performed in the bottle, with a minimum of buffer that can be easily disposed of. Buffer capsules have a long shelf life and are a convenient way of storing larger quantities. Professional and low ionic strength buffers are also available on request.

Buffers, Solutions and Accessories

Cat No.	Description
RLS010/MOD	Replacement pH sensing module
BS04/500	pH 4 buffer, twin necked bottle. 500 ml
BS07/500	pH 7 buffer, twin necked bottle. 500 ml
BS10/500	pH 10 buffer, twin necked bottle. 500 ml
BC54	pH 4 buffer capsules, five vials of ten capsules
BC57	pH 7 buffer capsules, five vials of ten capsules
BC510	pH 10 buffer capsules, five vials of ten capsules
BC5RM	Mixed box buffer capsules, 20 x pH 4, 20 x pH 7, 10 x pH 10
RCES	Cantilever electrode stand
EMK001	Electrode maintenance kit for use with electrodes with silver /siler chloride reference
EMK002	Electrode maintenance kit for use with electrodes with calomel reference
S00V/250	Electrode storage solution 250 ml
ESB	Electrode storage bottles (for 12 mm diameter types), pack of 3

