



Pipette Stands

6 place and 3 place stands for Finnpiette Digital single channel pipettes, which consists of rugged molded supports and an aluminum crossbar.



Focus Pipette Stand

Focus Pipette stand with clear polycarbonate cover, holds 6 pipettes



Novus Pipette Stand

Novus mechanical stand is designed to fit both single channel and multichannel Novus pipettes. The stand can also fit all other Finnpiettes as well eg Finnpiette Focus Plus Multichannels.



Multichannel Pipette Stand

The dual purpose Finnpiette multichannel stand provides convenient storage for 4, 8, 12 channel, digital, colour or varichannel multichannels.



Convenient Shelf Hanger

Every Finnpiette Digital and Focus comes with a convenient stick-on shelf hanger. So you can keep your pipette wherever you need it.



Safety Label

For easy identification and added security, all Finnpiette Digitals have a new safety label on top of the grippy finger rest. The clear plastic window covers an adhesive label which you can mark with your name or any other vital information.

Accessories for Multichannel Pipettes

Choice of three different reagent reservoirs:

- Standard autoclavable reagent trough (9510027).
- Base dimensions approximately 140 x 48 mm. Approximate capacity 60 ml. V-shape for maximum liquid recovery.
- 8 channel autoclavable reservoir (9510037).

- Standard microplate footprint. Capacity approximately 10 ml per channel. Sloped bottom for maximum liquid recovery.
- Autoclavable reservoir (9510047). As 9510037 above but with no channel divisions. Capacity approximately 125 ml.

