

BioMate has been especially designed for routine tasks in research and clinical laboratories where comfort and convenience are greatly appreciated. BioMate operates on the reliable ball-valve principle. BioMate can be used with all serological pipette brands covering a volume range of 1 ml to 100 ml.

BioMate's precision is ensured by the blow-out function that guarantees all liquid dispensed.

Product Features

- * Extremely easy and exact pipetting
- * Unique ergonomic design with thumb flow control
- * Special bench stand for storage and battery recharge
- * Lightweight for low stress injury risk

Specifications

Batteries	2 x 2/3 Aa NiMH
Battery Life Service	Approx - 8 hours
Charging Time	14 hours full discharged battery
Pipette Adaptor	Silicone
Filter	PTFE
Nose Cone Material	Polypropylene
Housing and stand	ABS
Lever and Coloured Ring	Polypropylene
Charger	2.8V 60 mA
Pipette Types	Plastic or Glass up to 100 ml, Pasteur Pipettes



Ordering Information

Cat No.	Description	Qty
4580530	BioMate	1

Filters and Accessories

Filters

4580560	Filter, non-sterile	5
4580570	Filter, sterile	5

Accessories

4580580	Silicone pipette adaptor	1
4580590	Base with base plate	1
4580600	Base plate	1
4580610	Nose cone	1