

The improved design of Finnpiette Focus Multichannel Plus, combined with super blow-out, ensures improved accuracy and precision. Finnpiette Focus Multichannel Plus pipettes have the same ergonomic advantages as Finnpiette Focus single channel pipettes. The grippy finger rest lets the pipette balance at the optimum pipetting angle, improving pipetting results. Finnpiette Focus Multichannel Plus is much more comfortable and less tiring over long periods of repetitive pipetting. The soft-touch tip ejection is more efficient and requires less effort, which is especially important with multichannel pipettes.

Product Features

- 8-, 12- & 16-channel multichannel pipettes for faster, easier repeat dispensing
- Superior accuracy and precision
- Fast volume setting
- Ergonomic multichannel pipetting
- Soft-touch tip ejection
- Light weight

Applications

- Clinical
- Biotechnology
- Chemical
- Environmental
- Food



Ordering Information and Specifications

Cat No.	Channels	Volume Range	Increments	Accuracy		At	Precision s.d. * μ l	CV %	Colour Code	Finntip
				μ l	%					
4610100	8	1-10 μ l	0.1 μ l	\pm 0.24 \pm 0.12	\pm 2.4 \pm 1.2	10 μ l 1 μ l	0.16 0.08	1.6 8	pink	10
4610120	8	5-50 μ l	0.1 μ l	\pm 0.75 \pm 0.25	\pm 1.5 \pm 5	50 μ l 5 μ l	0.35 0.1	0.7 2	yellow	250 Univ 200 Ext
4610140	8	30-300 μ l	1 μ l	\pm 1.5	\pm 5.0	300 μ l 30 μ l	0.6	2	orange	300
4610110	12	1-10 μ l	0.1 μ l	\pm 0.24 \pm 0.12	\pm 2.4 \pm 1.2	10 μ l 1 μ l	0.16 0.08	1.6 8	pink	10
4610130	12	5-50 μ l	0.1 μ l	\pm 0.75 \pm 0.25	\pm 1.5 \pm 5	50 μ l 5 μ l	0.35 0.1	0.7 2	yellow	250 Univ 200 Ext
4610150	12	30-300 μ l	1 μ l	\pm 1.5	\pm 5.0	300 μ l 30 μ l	0.6	2	orange	300
4610170	16	5-50 μ l	0.1 μ l	\pm 0.75 \pm 0.25	\pm 1.5 \pm 5.0	50 μ l 5 μ l	0.35 0.10	0.7 2.0	turquoise	50

*s.d. = Standard Deviation CV = Coefficient of Variation