

Pipettes : Finnpiettes

Digital - Fully Autoclavable

"Snag-Free" Thumb Button Volume Setting

No danger of gloves becoming caught, positive click volume setting ensures no "wandering".

Unique Safety Label And Grippy Finger Rest

Ergonomic Body Design

Polymer construction gives an extremely robust pipette with low weight, good chemical resistance.

Individual Serial Numbers

Clearly marked on handle for traceability.

Unique "Super Blow-Out" Mechanism

A unique telescopic piston boosts the power of the blow-out, preventing capillary action and achieving efficient sample delivery, guaranteeing you precise delivery of micro volumes (0.2 - 50 µl)

Colour Coded Thumb Button

Easy Access Calibration Ring

Straight forward user-calibration using the tools provided.

Unique "Soft-touch" Tip Ejector

Ensures easy de-tipping every time.

Smooth Piston Action

Damped action means no risk of vortex contamination with Finntips.

Chemically durable materials

Calibration

All Finnpiettes are factory calibrated and certified at minimum and maximum volumes

