



Refrigerators and Freezers

Upright Freezers Temprow

Upright Freezer Temprow S

Temperature range from -18°C to -28°C, 8 drawers with individual freezer compartments.

Features

Epoxy covered metallic external case.
Strong, molded, chemically resistant ABS interior.
Manual defrost with water evaporation system.
Door can be hung to open left or right.

Control panel

Mains power ON/ OFF switch.
Temperature thermostat regulator.
Temperature indicator lamp.
Audible and visual alarms.
Push button alarm reset/stop.
Push button with quick freeze indicator.



Upright Freezer Temprow M

Temperature range from -9°C to -26°C, 7 individual freezer compartments with 14 baskets.

Features

Epoxy covered external case.
Strong, chemically resistant, ABS molded case.
Specially coated 'Rilsan' wire mesh baskets.
Manual defrost, with water evaporation system.
The door can be hung to open left or right.
Lockable door.

Control panel

Temperature thermostat regulator.
Digital temperature display.



Upright Freezer Temprow L

With fan assisted air circulation. Temperature range from -15°C to -35°C, comes with 4 shelves

Features

External and internal case made completely of AISI 310 stainless steel, with water pump and drain.
Chemically resistant.
Acclimatised hermetically sealed CFC free compressor.
Automatic de-frost with electrical water evaporation system.
Shelves and rack.
Adjustable front level supports.
The door can be hung to open , left or right.
Door lock.

Control panel

Mains power ON/ OFF switch.
Temperature thermostat regulator.
Digital temperature display.



Part No.	Capacity litres	Range °C	H x W x D (internal) cm	H x W x D (external) cm	No. of drawers	Power W	Weight Kg	
C/R/020027	300	-18 to -28	158 46 43	184 60 63	8	110	74	Comes with drawers
C/R/020028	520	-9 to -26	150 64 53	172 75 72	14	345	111	Comes with baskets
C/R/020029	650	-15 to -35	144 62 67	199 76 70	4	600	160	Comes with shelves