

Refrigerators and Freezers

Chest Freezer Combatemp

Temperature range from -10°C to -30°C or -40°C



Features

Epoxy covered external case. Interior of aluminum, with guide supports for baskets or trays. The interior of the door is made of easy to clean, strong, chemically resistant, molded ABS.

Drop down lockable door with compensating support hinges. Silent, vibration and CFC free hermetically sealed compressor.

Automatic interior light when the door is opened.

Digital electronic temperature control. N.T.C. temperature probe.

Control Panel

Temperature control.

Digital live and preset temperature display.

Push-button temperature selector.

Part No.	Capacity litres	Temperature until °C	H x W x D (internal) cm	H x W x D (external) cm	Power W	Weight Kg
C/F/020001	250	-30	65 x 90 x 45	90 x 113 x 74	160	65
C/F/020002	352	-30	65 x 113 x 50	90 x 137 x 79	200	75
C/F/020003	250	-40	65 x 90 x 45	90 x 113 x 74	230	70
C/F/020004	352	-40	65 x 113 x 50	90 x 137 x 79	400	80