

LC3 Winchester Cabinet

The free-standing Winchester Cabinet has been designed to fit conveniently under the standard laboratory bench and will hold eight 80oz Winchester bottles.

Constructed from heavy-gauge sheet steel, the cabinet is attractively finished in stoved enamel and clearly marked with the International Standard Fire Warning Symbol.

Capacity: Two x four, 80oz Winchester bottles



LC3 Winchester Cabinet



DC1 Cabinet

DC1 Cabinet

This carefully designed cabinet incorporates several features to ensure secure storage of scheduled drugs and poisons.

The cabinet is intended for wall mounting direct through the specially reinforced panels incorporated in the back, and is fitted with a six-level security lock.

The three shelves may be adjusted at 75 mm centers.

Ordering Information

Cat No.	Dimensions (H x W x D) mm		Weight kg	No. of Shelves
	External	Internal		
LC3	705 x 350 x 310	680 x 320 x 300	18	1
DC1	610 x 455 x 155	605 x 450 x 150	16	3

Medical Cabinets

This heavy duty medical cabinet is ideal for storage of drugs and medical supplies in institutional and commercial situations.

Product Features

- Fixing holes in back only
- Adjustable removable shelves
- Smart easy clean finish

Ordering Information

Cat No.	Dimensions (H x W x D) mm		Weight kg	No. of Shelves
	External	Internal		
DC2	610 x 305 x 229	560 x 255 x 220	20	3
DC3	610 x 610 x 229	560 x 560 x 220	35	3