



Digital Vacuum Gauge



Deluxe Convenience Cart



Vacuum Pump Lubricating Oil



CryoCool Heat Transfer Fluid



Vacuum Pump Oil Filter

### Ordering Information

Cat No.	Description
DVG50-230	Digital vacuum gauge, complete with pressure transducer, 230 V/50/60 Hz
DVG50-115	Digital vacuum gauge, complete with pressure transducer, 115 V/50/60 Hz
CC120/DX	Deluxe convenience cart - 61 x 91.4 x 96.5 mm
SP01-B	One litre of vacuum pump oil
SP01	Case of 12 one litre bottles of vacuum pump oil
SP01-SB	Vacuum pump lubricating oil, synthetic, 1 litre
SP01-S	Vacuum pump lubricating oil, synthetic, 1 litre, case of 12
SFF1-B	One litre of vacuum pump flushing fluid
SFF1	Case of 12 one litre bottles of vacuum pump flushing fluid
VPOF110	Vacuum pump oil filter for models VLP80, VLP120, VLP200 and VLP285
A0310	Replacement oil filter cartridge for VPOF110
MF190	Mist filter replacement cartridge for OME190 oil mist eliminator

### Chemical Trap Kits

A chemical trap is highly recommended to absorb residual volatile vapours that may bypass the primary cold trap. Chemical trap kits are designed to efficiently and easily trap vapours and reduce emissions. Replaceable cartridges install easily inside the trap for quick changes once the trap becomes saturated. Dry chemical traps (SCT120 and cartridge) are installed between the cold trap and the oil-filled high-vacuum pump to help prevent pump damage from contaminating solvent vapours. Solution-based traps (ANT100 plus solution) and dry chemical traps are used on the exhaust of oil-free pumps to reduce emissions.



Chemical trap kit DTK120R

### Ordering Information

Cat No.	Description
DTK120R	Chemical trap kit, for absorbing volatile radioactives, includes activated carbon cartridge and all the fittings
SCT120	Chemical trap (order DC120A or DC120R cartridges separately)
DC120A	Disposable cartridge, with colour indicator for trapping acids and water vapour
DC120A/4	Disposable cartridge, with colour indicator, for trapping acids and water vapour (package of 4)
DC120R	Disposable activated charcoal cartridge, for trapping radioactivity
DC120R/4	Disposable activated charcoal cartridge, for trapping radioactivity (package of 4)
ANT100	Ammonia neutralising trap assembly for SPD1010, SPD2010, ISS110, UVS400A and UVS800DDA
ANS121/4	Ammonia neutralising solution (4 bottles) for use with ANT100
NSA300TF/4	Acid neutralising solution (4 bottles) for use with ANT400
NSA300TF	Acid neutralising solution (1 bottle) for use with ANT400